

WHEELCHAIR BRACKET



Dock N Lock is compatible with many wheelchair models as listed below. Unwin in conjunction with Elap is working continuously to increase the range so please contact Unwin Sales office, should your wheelchair not be listed.

MAKE	MODEL
Blank Plough Bar	DNLB10
Long Blank Plough Bar	DNLB35
DAYS HEALTHCARE	Escape PWR SE (DNLB9)
	Viper (DNLB33)
	Viper Plus (DNLB33)
EASYCARE	Genie (DNLB14)
HANDICARE	Alex (DNLB24)
	Beatle (DNLB11)
	Ibis Forward (DNLB25)
	Puma (DNLB11)
INVACARE	Harrier Plus (DNLB12)
	Mirage (DNLB6)
	Phoenix (DNLB13)
	Phoenix Plus (DNLB13)
	Pronto 61 (DNLB3)
	Spectra Blitz (DNLB18)
	Spectra Plus (DNLB8)
	Storm 3 (DNLB4)
	Storm 4 (DNLB36)
	TDX SP (DNLB2)
	Typhoon 2 (DNLB7)
XTR (DNLB23)	

MAKE	MODEL
LEVO	Combi Stand up (DNLB38)
KARMA	KP25 (DNLB17)
PRIDE MOBILITY	Quantum 600 (DNLB1)
	Quantum 600XL (DNLB1)
	Quantum 6000 (DNLB22)
Transport Unoccupied	Jazzy 600 (DNLB1)
Transport Unoccupied	Jazzy 600XL (DNLB1)
Transport Unoccupied	Jazzy Ultra 1103 (DNLB21)
SUNRISE	Quickie F55 (DNLB16)
	Groove FWD (DNLB26)
	Groove RWD (DNLB28)
	Salsa RWD(DNLB19)
	Salsa MWD(DNLB20)
	Rumba (DNLB5)

Note:

Pride Mobility: Please check with manufacturer to ensure that your wheelchair model configuration is transportable.

Sunrise Medical: wheelchairs fitted with Balle Seating Modules are not compatible with Dock N Lock.